



A Dräger and Siemens Company

Field Service Procedure

Part Number: SP00295

Rev: A

Date: 2 February 2004

© 2004 Dräger Medical

Innovian System Checkout Procedure

TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
1 PURPOSE	2
2 SOFTWARE REQUIREMENTS	2
3 RELATED DOCUMENTS	2
4 ACRONYMS	2
5 DEFINITIONS.....	3
6 PREREQUISITES.....	3
7 PROCEDURES.....	4
7.1 Server Test Procedures	4
7.1.1 Document General Information.....	4
7.1.2 Network Information.....	4
7.1.3 Database Information.....	5
7.1.4 Services Information.....	5
7.1.5 PCAnyWhere Verification	5
7.2 Workstation Test Procedure	7
7.2.1 Document General Workstation Information	7
7.2.2 Testing the Innovian System	8
7.2.3 Modifying Drug Entries.....	9
7.2.4 PCAnyWhere Verification	9
8 RESULTS SHEETS.....	10
8.1 SERVER RESULTS SHEET.....	10
A copy of this report shall be forwarded to Regulatory Affairs.....	10
WORKSTATION RESULTS SHEET.....	11
A copy of this report shall be forwarded to Regulatory Affairs.....	11

1 PURPOSE

This set of procedures provides instructions to test an Innovian Server or Workstation installation. The operator must have a basic understanding of Windows.

This series of procedures are broken into different sections for the (1) Server and (2) Clinical and Nonclinical Workstation.

2 SOFTWARE REQUIREMENTS

N/A

3 RELATED DOCUMENTS

<i>DOC NUM</i>	<i>DESCRIPTION</i>
Innovian-IP-IN-0322 (SP00292)	Innovian Server Installation Procedure
Innovian-IP-IN-xxxx (SP000293)	Innovian Workstation Installation Procedure (Upgrade, Install, and Disaster Recovery)

4 ACRONYMS

<i>ACRONYMS</i>	<i>DESCRIPTION</i>
DD	Designated Drive
ISU	Installation and Setup Utility Script
LAN	Local Area Network
NAD	North American Dräger
WPU	Workstation Processing Unit

5 DEFINITIONS

<i>DEFINITION</i>	<i>DESCRIPTION</i>
Click	When using a mouse, press the left mouse button. When using a touch screen, touch the item on the screen.
Designated Drive (DD)	The drive where the install is being done. Typically it is the C: drive unless designated otherwise .
Double-click	When using a mouse, rapidly press the left mouse button twice. When using a touch screen, rapidly touch the screen twice.
File	The contents of a folder.
Folder Tree/Directory Tree	The graphic representation of folders and subfolders.
Panel	When a screen is divided into two sections (left and right).
Right-click	When using a mouse, press the right mouse button. When using a touch screen, touch the item on the screen. After the item is selected, press and hold the Shift key, and then press F10.
Subfolder	A folder within a folder.

6 PREREQUISITES

Before this procedure is performed, the operator must obtain or access the following information:

- Innovian Workstation Installation Procedure (Upgrade, Install, and Disaster Recovery)
- Innovian Server Installation Procedure.

7 PROCEDURES

7.1 Server Test Procedures

7.1.1 Document General Information

1. Record the **Hospital Name**, **Address**, **Install Date** and **Server Location** on the Server Results Sheet.
2. Record the **Serial Number** and **NAD Asset Tag Number (if available)** on the Results Sheet.
3. Record the **Domain Username** and **Password** on the Results Sheet.
4. Record the **Local Administrator Username** and **Password** on the Results Sheet.

7.1.2 Network Information

1. Click **Start -> Run**.
2. Type **cmd** in the edit box and press the **Enter** key. The Command window appears.
3. At the Command prompt, type **cd ** and press the **Enter** key.
4. At the Command prompt, type **ipconfig /all** and press the **Enter** key. The IpConfig window appears.

```

C:\WINNT\System32\cmd.exe
Microsoft(R) Windows NT(TM)
(C) Copyright 1985-1996 Microsoft Corp.
C:\>ipconfig /all

Windows NT IP Configuration

    Host Name . . . . . : 
    DNS Servers . . . . . : 152.3.78.1
                          152.3.78.2
    Node Type . . . . . : Hybrid
    NetBIOS Scope ID. . . . . : 
    IP Routing Enabled. . . . . : No
    WINS Proxy Enabled. . . . . : No
    NetBIOS Resolution Uses DNS : No

Ethernet adapter NI002:

    Description . . . . . : Compaq Ethernet/FastEthernet or Gigabit NIC
    Physical Address. . . . . : 00-02-A5-6B-84-EC
    DHCP Enabled. . . . . : No
    IP Address. . . . . : 152.3.157.122
    Subnet Mask . . . . . : 255.255.0.0
    Default Gateway . . . . . : 152.3.157.250
    Primary WINS Server . . . . . : 152.3.173.246
    Secondary WINS Server . . . . . : 152.3.173.247

Ethernet adapter NI001:

    Description . . . . . : Gigabit Ethernet Adapter from Compaq
    Physical Address. . . . . : 00-40-A5-B1-A8-0E
    DHCP Enabled. . . . . : No
    IP Address. . . . . : 192.168.18.1
    Subnet Mask . . . . . : 255.255.255.0
    Default Gateway . . . . . : 
  
```

Figure 1 Server Command Window

5. Record the following information on the Server Results Sheet:

• <i>Host Name as Computer Name</i>	• <i>First DNS Server as DNS1</i>
• <i>Physical Address as NIC Address</i>	• <i>Second DNS Server as DNS2</i>
• <i>IP Address</i>	• <i>Primary WINS Server as WINS1</i>
• <i>Subnet Mask</i>	• <i>Secondary WINS Server as WINS2</i>
• <i>Default Gateway</i>	
6. At the Command prompt, type **explorer** and press the **Enter** key. The Windows Explorer window appears.

7. Right-click on **My Computer** and choose **Properties**. The System properties dialog appears.
8. Record the **Windows** version on the Results Sheet.
9. Select the **Network Identification** tab.
10. Record the **Domain** name on the Results Sheet.
11. Close the System properties dialog.
12. Close Windows Explorer and the Command window.

7.1.3 Database Information

1. Click **Start > Programs > Accessories > Explorer**.
2. Browse to **D:\Program Files\Microsoft SQL Server\MSSQL\Binn\sqlservr.exe**.
3. Right-click on **sqlservr.exe**.
4. Choose **Properties**.
5. Select the **Version** tab.
6. Record the **File Version** on the Results Sheet.
7. Select **OK** on the **Properties** dialog and close Explorer.
8. Click **Start > Programs > Microsoft SQL Server > Enterprise Manager**.
9. Browse to **Console Root > Microsoft SQL Servers > SQL Server Group > Local > Databases**.
10. Confirm the presence of: **Audit4sys, Patientinfo, Periopdm** and **Terminology** (Site Specific) databases.
11. Record the database names on the Results Sheet.
12. Expand the **Periopdm** database in Enterprise Manager.
13. Go to **Tables** and scroll down to the **VersionInfo** table.
14. Right-click on the version info table and select 'open table' and then 'return all rows'.
15. Record the latest entry in the **VersionInfo** table on the Results Sheet.
16. Close **Enterprise Manager**.

7.1.4 Services Information

1. Click **Start > Settings > Control Panel > Administrative Tools > Services**.
2. Confirm that the **Distributed Transaction Coordinator** service status is set to **Started**, and start-up type is set to **Automatic**.
3. Confirm **MSSQL Server** service status is set to **Started**, and start-up type is set to **Automatic**.
4. Confirm **SQL Server Agent** service status is set to **Started**, and start-up type is set to **Automatic**.
5. Successful start-up of these services will confirm that the databases are running.
6. Close the **Services** dialog.

7.1.5 PCAnywhere Verification

1. Click **Start > Programs > PCAnywhere** to open PCAnywhere.
2. Click the **Hosts** button to view all host items.
3. Double-click the host icon with the server name to start the PCAnywhere Host service.
4. Contact Saturn Support at (800) 523-6817 ext. 2362 to test the host connection.
5. If the **PCAnywhere connection is functional**, record the results on the Results Sheet.
6. Record **PCAnywhere Username** and **Password** on the Results Sheet.
7. On the Results Sheet, record the installed **version** of PCAnywhere and the type of installation as **H** for Host Only or **HR** for Host and Remote.

***** END OF THE Innovian SERVER TEST INSTRUCTIONS *****

7.2 Workstation Test Procedure

7.2.1 Document General Workstation Information

1. Record **Hospital Name**, **Address** and **Install Date** on the Workstation Results Sheet.
2. Record **Location** of the EG unit or WPU on the Results Sheet.
3. Record the anesthesia machine (to which Innovian is mounted) as **Mount** on the Results Sheet.
4. Record the type of **UPS** (WPUs only) it is connected to on the Results Sheet.
5. Record the type of workstation as **C** for clinical or **NC** for nonclinical on the Results Sheet.
6. For clinical workstations, record the **Type** of clinical as **EG** or **WPU** for Dräger or **HPU** for third-party clinicals.
7. Record the EG or WPU **serial number** and **NAD Asset Tag number** on the Results Sheet.
8. Record **Domain Username** and **Password** on the Results sheet.
9. Record **Local Administrator Username** and **Password** on the Results Sheet.
10. From the **Start** menu, click **Run**, type **cmd** in the edit box, and then press the **Enter** key.
11. At the command prompt, type **cd ** and press the **Enter** key.
12. At the command prompt, type **ipconfig /all** and press the **Enter** key. The IPConfig window appears.
13. Record the following information on the Workstation Results Sheet:
 - *Host Name as Computer Name*
 - *Physical Address as NIC Address*
 - *IP Address*
 - *Subnet Mask*
 - *Default Gateway*




Figure 2 Workstation Command Window

14. At the Command prompt, type **explorer** and press the **Enter** key. The Windows Explorer window appears.
15. Right-click **My Computer** and choose **Properties**. The **System Properties** dialog box appears.
16. Click the **General** tab.
17. Record the **Windows Version** on the Results Sheet.

18. Select the **Network Identification** tab.
19. Record the **Domain** on the Results Sheet.
20. Close the **System Properties** dialog by clicking **OK**.
21. Browse to **C:\Program Files\Microsoft SQL Server\MSSQL\Binn\sqlservr.exe**.
22. Right-click **sqlservr.exe** and choose **Properties**.
23. Select the **Version** tab.
24. Record the **File Version** on the Results Sheet.
25. Close the **Properties** dialog and **Explorer**.
26. From the **Start** menu, click **Settings** then **Printers**.
27. Record the name of each installed printer on the Results Sheet.
28. Right-click the printer icon and choose **Properties**. Record the **Driver** installed for that printer (displayed on the **General** tab).
29. Choose the **Ports** tab and record the port it uses. If LPR printing is utilized, record the IP address of the printer as well. (Repeat steps 27 and 28 for each installed printer.)
30. Restart the computer.

7.2.2 Testing the Innovian System

1. Does the System automatically log on (Site Specific) to Windows? Record the results on the Results Sheet.
2. Click **Start >Programs >Innovian Applications >Recorder**.
3. The Innovian Logon screen should appear. Enter the **User Name** and **Password**, then click **OK**.
4. Is the **Innovian logon successful**? Record the results on the Results Sheet.
5. Click **Help>About** on the Innovian screen and record the Innovian version information on the Results Sheet.
6. Click the **Log Out** toolbar button at the top of the screen. The Innovian Logon screen should appear.
(Steps 9 through 12 should be performed from one non-clinical or clinical workstation only)
7. At the 4ys Logon screen, enter the **User Name** and **Password**, then click **OK**.
8. Select the **Utilities** menu and choose **List Manager**. The List Manager window appears.
9. Click the plus sign (+) next to **Drugs** and select **Induction** in the left pane.
10. Deselect the check box next to the drug **thiopental** from the right pane.
11. Select the **File** menu and choose **Exit** to close List Manager. Save your changes.
12. Click the **Log Out** toolbar button at the top of the screen. The Innovian Logon screen should appear.
(Steps 16 through 24 should be performed from all WPU's, EGs, and nonclinical workstations)
13. At the Innovian Logon screen, enter the **User Name** and **Password**, then click **OK**.
14. In the Innovian Main window, select the **File** menu, and then choose **New Case**.
15. Click the **Drug** toolbar button at the top of the screen. The Add Drug dialog appears.
16. Click the **Induction** tab in the Add Drug dialog. Verify **thiopental** is not listed and click **OK**.
Note: If the above is true, **List Manager is functional**. Record the results on the Results Sheet.
17. Click the **Cancel** button in the Add Drug dialog box.
18. Click the **Close** toolbar button at the top of the screen.
19. A dialog box appears, asking if you want to save the case. Select **No**.
20. Click the **Log Out** toolbar button at the top of the screen. The Innovian Logon screen should appear.

7.2.3 Modifying Drug Entries

(Perform this step from one WPU, EG or non-clinical workstation only)

1. At the Innovian Logon screen, enter the **User Name** and **Password**, then click **OK**.
2. Select the **Utilities** menu and choose **List Manager**. The List Manager window appears.
3. Click the plus sign (+) next to **Drugs** and select **Induction** in the left pane.
4. Select the check box next to the drug **thiopental** in the right pane.
5. Select the **File** menu and choose **Exit** to close List Manager. Save your changes.

7.2.4 PCAnywhere Verification

1. Click **Start > Programs > PCAnywhere** to open PCAnywhere.
2. Click the **Hosts** button to view all host items.
3. To start the PCAnywhere Host service, double-click the host icon with the workstation name.
4. If the **Host Service is functional**, record the results on the Results Sheet.
5. Record **PCAnywhere Username** and **Password** on the Results Sheet.
6. On the Results Sheet, record the installed **version** of PCAnywhere and type of installation as **H** for Host Only or **HR** for Host and Remote.
7. Restart the computer.

***** END OF THE Innovian WORKSTATION TEST INSTRUCTIONS *****

8 RESULTS SHEETS

8.1 SERVER RESULTS SHEET

Sales Order Number: _____

Hospital Name: _____ Install Date: _____

Street Address: _____

City: _____ State: _____ Zip: _____

Serial Number: _____	NAD Asset Tag: _____
NIC Address: _____	IP Address: _____
Computer Name: _____	Subnet Mask: _____
Domain: _____	Gateway: _____
Domain Username: _____	DNS1: _____
Domain Password: _____	DNS2: _____
Local Admin Username: _____	WINS1: _____
Local Admin Password: _____	WINS2: _____
PCAnywhere Username: _____	Server Database Name: _____
PCAnywhere Password: _____	Server Database Version: _____
	Server Location: _____

VERSION INFORMATION:

Windows: _____ SQL: _____ PCAnywhere: _____ PCA Type: _____

VERIFICATIONS:

Services are functional

PCAnyWhere connection is functional

All Verifications Passed Successfully Initials _____

NOTES:

Signature: _____ Date: _____

A copy of this report shall be forwarded to Regulatory Affairs

WORKSTATION RESULTS SHEET

Sales Order Number: _____
Hospital Name: _____ Install Date: _____
Street Address: _____
City: _____ State: _____ Zip: _____

Serial Number: _____ NAD Asset Tag: _____
NIC Address: _____ IP Address: _____
Computer Name: _____ Subnet Mask: _____
Domain: _____ Gateway: _____
Domain Username: _____ Server Database Name: _____
Domain Password: _____ DB User ID: _____
Local Admin Username: _____ Flat Panel Serial Number: _____
Local Admin Password: _____ Location: _____ PCA Type: _____
PCAnywhere Username: _____ Mount: _____ UPS: _____
PCAnywhere Password: _____ C/NC: _____ Type: _____

VERSION INFORMATION:

Windows : _____ SQL: _____ Innovian: _____ PCAnywhere: _____

VERIFICATIONS:

Version Labels are correct
Innovian logon is successful
List manager is functional
Innovian Recorder is functional
PCAnyWhere Host Service is functional

All Verifications Passed Successfully Initials _____

PRINTERS: (Record each installed printer's name, path, IP address and driver; Specify if LPR printing is utilized)

NOTES:

Signature: _____ Date: _____

A copy of this report shall be forwarded to Regulatory Affairs



A Dräger and Siemens Company

DrägerService is a division of
Draeger Medical, Inc.
3122 Commerce Drive
Telford, PA 18969
Tel: (215) 721-5402
(800) 543-5047
Fax: (215) 721-5784
Web: www.draegermedical.com
Printed in the U.S.A.